Peeling Ultrasonic

Skin hygiene and regeneration system through ultrasonic micro-vibrations. **Peeling Ultrasonic** removes the superficial layers of the skin producing a deep cleansing of the pores and an improvement in the appearance of the spots.

Peeling Ultrasonic is a technique to be included in facial and body beauty treatments that aim for deep hygiene, peeling or rejuvenation of the skin.

2 VERSIONS

PEELING DIAMOND STATION

Equipment with diamond tips for dermabrasion + ultrasonic peeling with ultrasonic micro vibrations for cleansing, skin regeneration, anti-marks and anti-wrinkle. Treatments contain 1 spatula for facial and body treatments + 6 rounded Diamond tips.

PEELING ULTRASONIC XTRA I

Ultrasonic micro vibration equipment for cleansing, skin regeneration, anti-marks and anti-wrinkle treatments

BENEFITS

- Exfoliates, moisturises and cleanses without eroding the skin, effectively and quickly.
- ✓ Stimulates the regeneration of new skin cells as it removes the superficial layers of the skin.
- Achieves a rejuvenated and luminous appearance from the very first session.
- ✓ The technique of Peeling Ultrasonic XTRA is not aggressive or invasive, it has no side effects to the skin.

Also: It naturally peels the skin, deeply cleanses the pore, rejuvenates the dermis, improves penetration of cosmetics and reduces the colour of skin spots.

Technical Characteristics

Power: 230-110 V / 50-60 Hz.

Consumption: 50 W

Weight: 5 Kg

Dimensions: 16.1 x 30.5 x 25 cm







Peeling Ultrasonic Xtra

COSMETICS



DESCALING TONIC 500ML.

Cleansing and descaling tonic lotion. Slightly astringent. Cod. A-8906



REGENERATING MIXTURE

Refreshing concentrated solution for the regeneration of all skin types.
Cod. A-8907-50



AZULENE DESCONGESTANT TONIC 500ML.

Azulene speciÿc tonic solution to prepare the skin for a complete treatment.

Cod. A-8909



SWEETYCARECREAM 500ML

Hydrating, regenerating and calming cream. indicated as a mask for post-peeling treatment. Peeling.

Cod. A8925-500