## UNIK DERMAJET



# DISCOVER OUR BIG NEW

#### THE NEW GENERATION OF TREATMENTS **OF BEAUTY TO TREAT BODY AND FACE**

UNIK is the exclusive multiplatform that works synergistically with different technologies focused on integral body remodeling and facial treatments, allowing the treatment of different skin conditions.

It is a device that is integrated into the range of successful Dermajet ® devices with a long history at the service of body remodeling. Now adding his technique in facials.

Unik models the figure, reduces the body contour and reaffirms the tissues, thanks to its technologies that work simultaneously, increasing its effectiveness.

#### **BODY TECHNOLOGIES**



With TRIFREQUENCY RADIOFREQUENCY we work at 1 Mhz - 740 Khz - 420 Khz cell stimulation at 3 different depths that allows the fibroblast to produce collagen (firmness), elastin (elasticity), hyaluronic acid (volume and hydration) as well as the elimination of localized fat. With the SEQUENTIAL **VACUUM** we work on a detachment massage that breaks cellulosic sclerosis by traction, increasing blood flow and opening the lymphatic channels with different wavelengths that will allow us to carry out a multitude of treatments. We work with ULTRASOUNDS at 1 MHZ that provide us with lipid capacity and sonophoresis. And we also work with LLLT - PHOTOBIODYNAMIC THERAPY with high power diodes.

Its multiple treatments are designed under our Morfologic System, thus managing to personalize the treatments.

#### **FACIAL TECHNOLOGIES**

Unik incorporates an exclusive triple modality handpiece for facial treatments and small areas.



- 🗸 Blood Drainage
- ✓ Sports
- Acne and rosacea
- Hyperpigmentation
- Scars

### LLLT PHOTOBIOTHERAPY

Low Lever Laser Therapy (LLLT) is a new system for treating different problems through light by localized light diodes, with certain wavelengths (colors) and certain working frequencies (Nogier).

#### I.O. 630 y 660 +- 10nm

- Cell regeneration
- Anti-inflammatory
- Elastin stimulation

- Stimulation of fibroblast and collagen production

#### I.O. 425 +-10 nm

- Bacteriological - Acts at the epidermal
- level

- Inhibits the growth of the horny layer

### COSMETICS

#### I.O. 530 +-10nm

- Inhibits excess melanin
- Anti-inflammatory
- Elastin stimulation
- Collagen stimulation

#### I.O. 830 +- 10 nm

- Relieve inflammation
- Treat insomnia
- Relieve pain, reduce
- stress and anxiety
- Activates circulation

E-HEBE

#### I.O. 590 +-10nm

- Rosacea
- -Dermatitis
- Anti-inflammatory
- -Lymphotropic
- Increases blood circulation
- Draining

#### NIX



Ultrasonic thermal and mechanical conductivity. They activate the elimination of fats, stimulate microcirculation and regenerate damaged tissues. E-HEBE

Professional cosmetics for rejuvenating and regenerating treatments. Formulated to stimulate the fibroblast by increasing the levels of collagen, elastin and hyaluronic acid.



Photodynamic booster Nix Wrinkles: Promotes neocollagenesis and significantly improves the damage caused by photoaging.

Nix Acne: for type I and II acne.

Nix Redness: Stimulates oxygenation and cellular respiration.

Nix Dark Spots: spots on the face and neck.



## Manufacturing beauty since 1985

C/ Frederic Soler, 57E 08224 - Terrassa (BCN) Teléfono: +34 93 184 43 73 enco@encoworld.com

www.encoworld.com/en